
- 1 Brush head cap
- 2 Brush head
- 3 Handle
- 4 Power on/off and mode button
- 5 Brushing modes
- 6 Post brushing feedback indicator
- 7 Battery indicator
- 8 Charger

Note: Power adapter not included. Charge only with IPX4 adapters with output ratings 5 V , ≥ 2 W. It is recommended to use a Philips adapter; see the Philips website for current information.

Locating the model number

Look on the bottom of the Philips Sonicare toothbrush handle for the model number HX52xx.

- 1 Align the brush head so the bristles point in the same direction as the front of the handle.
- 2 Press the brush head down onto the metal shaft until there is a small 1-2 mm gap.

Charging your Philips Sonicare

- 1 Place the power toothbrush handle on the charger

- 2 Connect the charger to a 5 V USB wall adapter.
- 3 Your toothbrush continues to charge until it is fully charged.
 - When the toothbrush is charging, the battery indicator pulsates green.
 - When fully charged, the battery indicator will light up solid green for 30 seconds before turning off.

Using your Philips Sonicare toothbrush

If you are using your Sonicare toothbrush for the first time, it is normal to feel more vibration than when using a non-electric toothbrush. It is common for first time users to apply too much pressure. Only apply gentle pressure and let the toothbrush do the brushing for you. Please follow brushing instructions below for the best experience.

If you are using your Sonicare toothbrush for the first time, it is normal to feel more vibration than when using a non-electric toothbrush. It is common for first time users to apply too much pressure. Only apply gentle pressure and let the toothbrush do the brushing for you. Please follow brushing instructions below for the best experience.

- 1 Wet the bristles and apply a small amount of toothpaste.
- 2 Place the toothbrush bristles against the teeth at a slight angle (45 degrees) (Fig. 7). The handle has two 45 - degree indication bars to help guide your Bass brushing technique. Press the brush head gently against the gum to make the bristles reach the gum line or slightly beneath the gum line.
- 3 Press the power on/off button to turn on the Philips Sonicare.

- 4 Gently keep the bristles placed on the teeth and in the gumline. Brush your teeth with small back and forth motion so the bristles reach between the teeth. Continue this motion throughout your brushing cycle.

- 5 To clean the inside surfaces of the front teeth, tilt the brush handle semi-upright and make several vertical overlapping brushing strokes on each tooth (Fig. 8).

6 Begin brushing in section 1 (outside upper teeth) and brush for 30 seconds before you move to section 2 (inside upper teeth). Continue brushing in section 3 (outside lower teeth) and brush for 30 seconds before you move to section 4 (inside lower teeth) (Fig. 9).

- After you have completed the brushing cycle, you can spend additional time brushing the chewing surfaces of your teeth and areas where staining occurs. You can also brush your tongue, with the toothbrush turned on or off, as you prefer (Fig. 10).

Battery status

This Philips Sonicare toothbrush is designed to provide at least 120 brushing sessions after fully charged for 10 hours, each session lasting 2 minutes. When the battery is low, the battery indicator starts to flash amber.

Note: To keep the battery fully charged at all times, you may keep your Philips Sonicare toothbrush on the charger when not in use.

You should clean your Philips Sonicare toothbrush regularly to remove toothpaste and other residue. Failure to clean your Philips Sonicare toothbrush may result in an unhygienic product and damage to it may occur.

Warning: Do not clean with boiling hot water.

If you are not going to use the product for a long time, unplug it, clean it (see chapter 'Cleaning') and store it in a cool and dry place away from direct sunlight.

Replace Philips Sonicare brush heads every 3 months to achieve optimal results. Use only Philips Sonicare replacement brush heads.

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